

# AMT™ Vial Crimp Sealing Machines

Pharmaceutical grade  
stainless steel construction

Up to 300 vials per minute  
production speeds

Stainless steel vibratory feeder

Tooling with quick disconnects  
for swift changeover

Minimal foot print

Crafted for lasting performance

Low maintenance



AMT™ Fully Automatic Crimp Sealing



Non destructive and objective  
crimp seal tester



Crimp sealing option in multi-  
functional AMT™ Servo Monoblock

Automated Machine Technologies, Inc.  
10404 Chapel Hill Road, Unit 100  
Morrisville, NC 27560-1186

Office@AMTLiquidFilling.com  
www.AMTLiquidFilling.com  
Phone: (919) 361 0121  
FAX: (919) 481 2120




AMT™ Semi-Automatic Crimp Sealing

- Standalone Tabletop Model
- Stainless steel, Teflon coated hard anodized aluminum and Delrin™ construction
- Non-shearing crimp sealing, clean operating; no crimping rollers are used
- Operator manually places crimp seal over the stoppered vial and pushes vial into the input side of the Crimp Sealing Machine
- Vial is automatically transferred through Crimp Sealing Machine via motor driven urethane belt and is simultaneously crimped by non-shearing blade; urethane belt minimizes particulate matter generation
- Crimped vials exit output side of Crimp Sealing Machine and accumulate on table

AMT™ has also developed a none destructive and objective method for testing the firmness of the crimped seal.

To see a video of various AMT™ models in production, please visit our website: [www.AMTLiquidFilling.com](http://www.AMTLiquidFilling.com)

## AMT™ Vial Crimp Sealing Machine Dimensions



Width=36", Depth=24", Height=52", Weight: 250 lbs

Power: 110V/6A or 250V/3A, 50-60Hz

As part of our policy of continuous product development we reserve the right to alter specifications without notice.

Automated Machine Technologies, Inc.  
10404 Chapel Hill Road, Unit 100  
Morrisville, NC 27560-1186

Office@AMTLiquidFilling.com  
www.AMTLiquidFilling.com  
Phone: (919) 361 0121  
FAX: (919) 481 2120